**CEREPLAS**

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---

# **CEREFORM<sup>®</sup>**

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---

# USER LEAFLET - CEREFORM®

## SILICONE BREAST IMPLANTS PRE-FILLED WITH SILICONE GEL

### I. DESCRIPTION OF THE DEVICE

CEREFORM® breast implants are implantable, long term, silicone elastomer medical devices pre-filled with silicone gel. These prostheses are sold sterile after ethylene oxide sterilisation. They are for one patient only, single use.

CEREFORM® breast implants are placed in a double packaging that ensures a dual microbiological barrier guaranteeing sterility until implantation. The cardboard box provides additional protection to the double packaging to deliver the product in the best possible conditions.

The traceability of each implant is assured by a unique identification number etched on its patch. The accompanying documents also carry this number and all the characteristics related to the implant.

CEREFORM® breast implants are manufactured with the latest generation of medical grade silicones which guarantee the necessary properties in order to be in conformity with the regulation requirements in force. An integrated barrier layer strongly limits silicone gel bleeding through the shell. Moreover, the specific formulation of the silicone gel makes the implant very similar to a natural breast. CEREFORM® breast implants do not contain any phthalates.

CEREFORM® breast implants exist in two shapes:

**The round shape** gives volume to the upper part of the breast and uplifts the bust. It is particularly well suited to women that already have formed breasts.

**The anatomical shape** gently follows the bust contour without a rupture between the upper breast and the thorax. It is ideal for very slim figures and women with small or no breasts.

Depending on their shape, CEREFORM® breast implants are available with two different surfaces:

**The smooth surface**, available in the round shape, for easy implantation and removal.

**The intermediate texture**, available in the round and anatomical shapes, for the best cell colonisation and a reduced risk of the appearance of retractile capsular contracture, whilst preserving ease of implantation and removal.

All these implants are available in different sizes. To select the best suiting reference, CEREPHAS® has developed a range of sizes identical to that of the implants. Refer to the literature for further information.

### II. ADVANTAGES

- High resistance to tearing and rupture thanks to the latest generation of medical silicone;
- Individual identification of the implant by laser etching on the patch allowing traceability for the whole implant lifetime;

- A strong-resisting patch, undetectable to the touch that integrates a silicone gel barrier layer to limit the effects of gel bleeding;
- Natural breast consistency thanks to the optimized gel silicone formulation;
- Strong, ergonomic double packaging designed for simple use and maximum protection;
- Two surfaces, which offer the choice between either easier implantation and explantation (smooth surface) or reduced capsular contracture prevalence (intermediate texture surface).

### **III. INDICATIONS**

CEREFORM® breast implants are indicated for:

- Mammary reconstruction following a mastectomy;
- Breast enlargement for aesthetic reasons;
- The correction of various acquired or congenital defects: dissymmetry, amastia, aplasia, hypomastia, hypoplasia;
- Replacement after explantation of an old or defective implant.

CEREFORM® breast implants are supplied sterile. They are not to be re-used or re-sterilized.

The practitioner is responsible for the pre-operative evaluation of the patient and the operating method used. He must inform his patient of the risks related to the intervention and eventual post-operative complications. The practitioner has to obtain the patient's informed consent by his signature on the specific document before the intervention. He also has to inform the patient of alternative methods to implantation.

### **IV. CONTRAINDICATIONS**

The implantation of CEREFORM® breast prostheses is subjected to the following contraindications:

- Pre-existing pathology in the site of implantation;
- Breast cancer in progression;
- General or local infectious condition;
- Pregnancy or breast-feeding in progress;
- History of autoimmune disease;
- Repeated unsuccessful implantation of the same kind of implants;
- Treatment by irradiation, microwave diathermy or steroids in progress;
- Known hyper sensibility to silicone;
- Psychological instability;
- Tissue insufficiency.

### **V. POSSIBLE COMPLICATIONS**

Like all surgical operations, the implantation of breast prostheses entails operative and post-operative risks. Every patient's body reacts differently to implantation. Before the operation, the patient should be informed of the eventual risks and complications related to the intervention and the introduction of the implant.

#### **Risks related to the surgery**

The introduction of mammary implants entails risks inherent to the surgery itself, such as the risks and complications of general anaesthesia. It is important to take into account all these risks in the pre-operative check-up and to inform the patient before the intervention.

### **Serous effusion without infection**

An effusion of serous fluid may occur because of the intervention or because of a traumatism in the region of implantation. It is necessary to watch over the after-effects of all traumatism in the region of implantation. Effusion can be treated with the help of a drain or a puncture while taking all necessary precautions to avoid damaging the implant. If effusion continues, removal of the implant can be foreseen.

### **Haematoma or oedema in the site of implantation**

A meticulous haemostasis should be done during the intervention to prevent haematoma in the site of implantation. In the event of persistent haematoma, a puncture can be done while taking all necessary precautions to avoid damaging the implant. Appropriate medical immobilization in the site of implantation during the weeks following the intervention allows reduction of post-operative oedema.

### **Loss of sensation in the site of implantation**

A partial or total loss of sensation in the nipple area is possible. These changes may be temporary (with a progressive return to normal in a year) or definitive. No treatment is available to date.

### **Post-operative pain**

Patient dependant intensity post-operative pain due to the surgical procedure is felt within two or three days following the intervention. Discomfort may continue during the following month. This pain can be treated with analgesics. Any persistent pain in the site of implantation should be examined in order to eliminate an eventual complication.

### **Asymmetry**

Post-operative asymmetry is the consequence of an incorrect choice of implant (size, shape), disproportionate in comparison with the other breast or a tissue reaction different from one breast to the other. If the asymmetry is important and causes the patient's discontent, removal or replacement of the implant can be foreseen. Asymmetry appearing several months or years after the implantation could mean capsular contracture or a rupture of the implant. In this case, a detailed examination is necessary and the prosthesis explantation can be foreseen.

### **Unaesthetic inflammatory scars**

Poor healing either natural or as a result of a complication may lead to the appearance of unaesthetic hypertrophic scars or keloids. Preventive or curative treatment of these scars using silicone dressings is advised. If the problem persists, surgical operation on the scars can be foreseen.

### **Post-operative infection**

Short or long term post-operative infection is very uncommon within the framework of mammary prostheses implantation. However, any infection should be treated when encountered. If the anti-biotherapy does not enable treatment of the infection, the implant can be removed.

### **Capsular contracture formation**

The scar envelope, which forms naturally around any foreign element implanted in the human body, can retract around it, compressing it abnormally. This painful retraction could lead to deformation of the breast and the rupture of the implant. Removal of the prosthesis (with or without re-implantation) can be foreseen. Squeezing is highly discouraged to treat this complication because of the risk of rupture.

### **Necrosis of adjacent tissues**

Tissue necrosis can be induced by:

- The implant soaking in an iodine solution before the implantation ;
- Abnormal local tissue reaction, which could be due, for example, to an infection, or to a radiotherapy tissue treatment performed before implantation, etc ;
- Significant tissue tension due to tissue insufficiency or to a too large implant.

Tissue necrosis can lead to poor cutaneous healing of the tissues next to the breast implant. The severity of the necrosis depends on its spread. If this necrosis has a functional impact or is painful, removal of the implant should be foreseen.

### **Rupture of the implant**

The implant shell may rupture following an operative traumatism (damaged implant during insertion or by surgical instruments), or a post-operative traumatism (violent shock, excessive compression of the mammary region) or due to implant natural ageing. This rupture could be asymptomatic or display deflation of the prosthesis or a change in the shape or appearance of the breast. In case of doubt, it is necessary to conduct a diagnostic examination (mammography, ultrasound or MRI) in order to ensure that the implant is in good condition. Lastly, regular follow-up enables early detection of an eventual rupture. If the rupture is confirmed, removal of the implant is necessary. To guarantee constant mechanical properties and limit the risk of rupture, tests in compliance with the actual standards are regularly performed.

### **Displacement of the implant**

Displacement of the implant can result from poor initial positioning, an injury in the area of the implantation or a significant and premature slackening of the surrounding tissue, which is insufficient to maintain the implant. This results in loss of implant functionality (herniation of the implant, implant turnaround or change in the shape of the breast), which requires a repeat intervention.

### **Deflation of the implant**

Deflation of the implant occurs if it is ruptured. This phenomenon is rare with silicone gel implants because of the cohesive nature of the gel. Any deflation perceived by the patient must be interpreted as a rupture of the implant and should lead to more in-depth examination. Any confirmed rupture of the implant requires its explantation.

### **Implant perceptible to the touch**

Poor initial positioning of the implant, inappropriate size, displacement of the implant or a too thick scar envelope around it, can make the implant perceptible to the touch. If this causes the patient to be unhappy or if retractile capsular contracture occurs in the capsule around the implant, a new surgery can be foreseen.

### **Ptosis of the breast**

Just like a natural breast, a breast with a mammary implant can develop ptosis over years due to tissue distension in the site of implantation. Ptosis is not dangerous. It can be surgically corrected.

### **Wrinkles / folds / protuberances / waves in the area of the implant**

It is possible that the shell of the implant creases or curls, forming waves depending on its holding in the pocket and the position of the implant in relation to the pectoral muscle, depending on the surgical indication. The creases can be perceptible on the surface of the skin. Only explantation can correct this phenomenon.

### **Calcification of tissues surrounding the implant**

Calcification is the calcareous deposit in the tissues adjacent to the mammary implant. These deposits are painful and may damage the prosthesis that will then be explanted. This is not a very widespread phenomenon.

### **Siliconomas, granulomas**

Siliconomas form when small quantities of silicone spread at a distance from the site of implantation. The body then creates a thin scar envelope around the silicone. Siliconomas are not dangerous but are indicative of a rupture or major gel bleeding of the implant. If a rupture is confirmed, the implant should be removed.

### **Perspiration / bleeding / leakage of silicone gel**

Despite its barrier effect, the silicone shell is not perfectly silicone gel leakproof. Therefore, small quantities of silicone may diffuse out of the implant and spread in the scar envelope as well as in the tissues. Silicone gel is not toxic for the organism but local reactions with formation of small scar capsules may occur.

### **Replacement of the implant / revision surgery**

The limited life of the implant and the above mentioned risks could lead to a revision surgery in order to maintain the desired effect. The patient should understand and accept the risks of revision surgeries before deciding on the implantation.

### **Systemic complications and autoimmune diseases**

Presently, there is no proof of correlation between the appearance of systemic complications, cancer or autoimmune diseases and silicone gel breast implants.

### **Definitive removal of the implant without replacement**

If various complications associated with the implant should recur, or if the surgeon judges that the patient's state of health requires the withdrawal of the implant, definitive removal without the possibility of replacement of the implant and the associated inaeesthetic result (sagging chest, wrinkles) must be considered.

## **VI. PATIENT WARNINGS**

Breast implants have a limited lifespan. Given the possible complications mentioned above and the natural wear and tear on the implant in the body (daily mechanical strain), the implant is likely to require withdrawal or replacement, which may involve another surgical intervention. The practitioner should warn their patients of all the complications that could occur during surgery and after the implantation and inform them of the alternatives to the implantation of breast prostheses (wearing of external prostheses, breast reconstruction, etc.).

The patient must be aware that if he has his implant(s) removed without re-implantation, the result will be unaesthetic (sagging breast, wrinkles, etc.). In order to proceed with the implantation, the practitioner should obtain the patient's free and informed consent using the specific document. Practitioners are advised to insist on the importance of normal checks for the detection of breast cancer, in addition to post-operative follow-up in order to detect any possible complications. The patient should comply with inspection visits and report any traumatism, deflation or pain in the area of the implant.

It is also essential to insist on the precautions the patient should take, particularly the need to keep with them their completed patient card provided by the practitioner, the use of topical drugs near the breasts, the need to point out the presence of a breast implant during imaging tests or medical interventions in the area close to the implantation, the need to wait at least three months after the intervention before getting pregnant and the importance of self-examination of the breasts every month. For the latter, we recommend showing the patient how to do this.

## **VII. METHOD AND FREQUENCY OF EVALUATION OF THE IMPLANTED PROSTHESIS**

A clinical examination of the patient and of the implant is recommended 1, 3, 6 and 12 months after the implantation of the prosthesis and then every year even if there are no symptoms related to the implant.

In case of deflation of the implant, traumatism sustained since the last examination, pain, deformation of the implant or any other sign that could be suspected capsular retractile contracture or a rupture of the implant, a complementary examination with imaging by mammography, ultrasonography or MRI should be done. A yearly examination with imaging for the detection of ruptures is also advised beginning with the eighth year of the implantation.

Since silicone is a partially radio-opaque material, the implant can obscure the zone subjacent to the implant during a medical imaging test (mammography, echography, MRI). The radiologist will need to use a different technique in order to choose the appropriate angles of view to compensate for potential overshadowing by the implant of the areas to be observed. For mammography, it is imperative that the technique is adapted, as there is also a risk of rupture or weakening of the implant during the examination.

According to medical literature based information, the estimated average lifetime of silicone gel filled breast implants is 12.9 years (bibliographical study made by CEREPHAS®). However, taking into account the activities and numerous parameters which may influence the life of the implant, it is not possible to give a precise estimate of the life of a breast implant for any given patient.

## VIII. USE OF THE IMPLANT

CEREFORM® breast implants must be used only by competent practitioners that are experienced in breast surgery. CEREFORM® breast implants are to be used in an operating room.

- Check that the double blister is in perfect condition and that the expiration date has not been reached. In case of doubt concerning one of the two wrappings, do not use the implant;
- Handle the implant under aseptic conditions;
- The implant should not be in contact with any foreign substance (talc from gloves, tissues, iodine solution, etc.);
- It is possible to immerse the prosthesis in a bath of sterile physiological saline solution brought to body temperature to facilitate the implantation;
- The integrity of the device should not be compromised in any way. The user must never add any liquid filling to the interior of the capsule.
- Be careful of sharp and pointed objects while handling the prosthesis;
- Check that the implant shell does not display creasing in its permanent position;
- The size of the incision should be adapted to the size of the implant to avoid a too large deformation of the implant which could deteriorate it and a wrench of the tissues due to important mechanical strains;
- Ensure proper positioning of the anatomical implant in its pocket using tactile orientation marks;
- Do not use medicated substances with the implant;
- The practitioner is advised to have an additional implant in stock in order to compensate for a defective implant or mishandling during the intervention;
- Do not massage or carry out any treatment by pricking (acupuncture) in the site of implantation after the surgery in order to avoid any deterioration of the implant;
- In case of unfortunate rupture of the shell, clean and rinse the implantation pocket with normal saline to eliminate all potential gel traces;
- Any damaged implant must be removed.

### Instructions for withdrawing the implant

- For withdrawing the implant, we advise removing the implant by using the same incision as that used for implantation;
- If it is found that the implant has ruptured, we strongly recommend cleaning the breast compartment with sterile physiological serum. If some of the gel from the implant is present in the compartment, it must be physically extracted by the surgeon.

### Disposal of the product

Any implant which has been removed or is not sterile must be disposed of with the biological waste.

## IX. TRANSPORT AND STORAGE OF THE IMPLANT

The CEREFORM® breast implant must be transported and stored in its cardboard wrapping, protected from exposure to light, at a temperature between 0 °C and 40 °C, and must be kept dry and in a horizontal position. It must be handled with care.

## X. RE-STERILISATION AND RE-USE

The CEREFORM® breast implant is supplied sterile for single use. The CEREFORM® breast implant is supplied in a sterile condition for single use. Any re-use of a CEREFORM® implant is categorically proscribed. Any re-use of a non-sterile CEREFORM® implant can cause

serious clinical complications which could be fatal. Furthermore, any resterilisation of a CEREFORM® implant is also categorically proscribed, as it may cause significant degradation of its mechanical properties and its dimensional characteristics.

## XI. LABELS AND PATIENT CARD

Each implant is supplied with six labels that group together information about the implant (manufacturer's details, reference and identification number of the implant). These labels should be completed with the patient's name, the practitioner's name and the date of the implantation and must be kept in the patient's medical records. A patient card is also supplied. A label will be attached to this card and the patient must keep it with her at all times to make it easier if there should be a need for any emergency medical care. One of the labels will also be stuck to the certificate signed by the surgeon, who attests of the patient's informed consent for implantation and its risks. This document is then returned to CEREPHAS® or to the authorized distributor in order to ensure complete traceability of the implant.

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